

BAUSCH & LOMB

**PLEASE READ CAREFULLY AND KEEP THIS PACKAGE INSERT FOR FUTURE USE IN CASE YOU
HAVE A PROBLEM. USE ONLY IF TAB AND SINGLE-USE CONTAINER ARE INTACT.**

ReNu® Preservative Free Lubricating & Rewetting Drops

For use with soft (hydrophilic) contact lenses, including disposable and lenses prescribed for frequent replacement.

DIRECTIONS:

Use as needed, 2-3 times a day or as directed by your eye care professional. Use while wearing lenses to moisten and lubricate your lenses and to relieve minor irritation, discomfort and burning. If discomfort persists **IMMEDIATELY** remove lens and contact your eye care professional.

- Detach one container from the strip by pulling tab toward you from top to bottom.
- Remove the liquid from the area around the opening by tapping.
- To open, twist the tab off the single-use container.
- During wear, put 2 or 3 drops on contact lens and blink several times.
- Discard container immediately after use. **DO NOT STORE AND REUSE.**

CONTENTS:

A sterile, borate buffered, isotonic solution containing benzadac lysine 0.25%, hydroxypropyl methylcellulose 0.25%, and potassium chloride in purified water.

ACTIONS:

ReNu Preservative Free Lubricating & Rewetting Drops, packaged in single-use containers, rehydrates, rewets, and helps to remove particulate material from the surface of the contact lens during lens wear and assists in maintaining normal lens hydration.

INDICATIONS (Uses):

ReNu Preservative Free Lubricating & Rewetting Drops, packaged in single-use containers, is indicated for use with soft (hydrophilic) contact lenses (including disposables and lenses prescribed for frequent replacement) as a lubricating and rewetting agent during the wearing period.

CONTRAINDICATIONS (Reasons not to use):

IF YOU ARE ALLERGIC TO ANY INGREDIENT IN THIS PRODUCT, DO NOT USE.

WARNINGS: PROBLEMS WITH CONTACT LENSES AND LENS CARE PRODUCTS COULD RESULT IN SERIOUS INJURY TO THE EYE.

It is essential that you follow your eye care professional's directions and all labeling instructions for proper use and care of your lenses and lens care products, including the lens case. Eye problems, including corneal ulcers, can develop rapidly and lead to loss of vision. Daily wear lenses are not indicated for overnight wear and should not be worn while sleeping. Clinical studies have shown the risk of serious adverse reactions is increased when these lenses are worn overnight. Extended wear lenses should be regularly removed for cleaning and disinfection or for disposal and replacement on the schedule prescribed by your eye care professional. Clinical studies have shown that there is an increased incidence of serious adverse reactions in extended wear contact lens users as compared to daily wear contact lens users. Studies have also shown that the risk of serious adverse reactions increases the longer extended wear lenses are worn before removal for cleaning and disinfection or for disposal and replacement.

Studies have also shown that smokers have a higher incidence of adverse reactions. If you experience eye discomfort, excessive tearing, vision changes, redness of the eye, immediately remove your lenses and promptly contact your eye care professional.

It is recommended that contact lens wearers see their eye care professional twice each year, or if directed, more frequently.

- To avoid contamination of your solution or your lenses:
Do not touch tip of container to any surface.
Do not store opened container, use immediately after opening and discard unused contents. **DO NOT STORE AND REUSE.**
- Not for use with heat (thermal) disinfection.

PRECAUTIONS:

- Always clean, disinfect, rinse and soak your lenses each time they are removed.
- Do not store opened container; use immediately after opening and discard unused contents.
- Store solution at room temperature.
- Use before the expiration date marked on the carton and container tab.
- Keep out of the reach of children.

ADVERSE REACTIONS (Problems and what to do):

The following problems may occur: eyes sting, burn or itch (irritation), comfort is less than when lens was first placed on the eye, feeling of something in the eye (foreign body, scratched area), excessive watering (tearing) of the eye, unusual eye secretions, redness of the eye, reduced sharpness of vision (poor visual acuity), blurred vision, rainbows or halos around objects, sensitivity to light (photophobia), or dry eyes.

If you notice any of the above:

- **Immediately remove your lenses.**
- If the discomfort or problem stops, then look closely at the lens.
- If the lens is in any way damaged, **DO NOT** put the lens back on your eye. Place the lens in the storage case and contact your eye care professional.
- If the lens has dirt, an eyelash, or other foreign body on it, or the problem stops and the lens appears undamaged, thoroughly clean, rinse and disinfect the lens, then reinsert it.
- If the problem continues, **IMMEDIATELY** remove the lens and consult your eye care professional.

If any of the above symptoms occur, a serious condition such as infection, corneal ulcer, neovascularization or iritis may be present. Seek immediate professional identification of the problem and prompt treatment to avoid serious eye damage.

All adverse reactions observed while using Bausch & Lomb ReNu Preservative Free Lubricating & Rewetting Drops should be reported to:
Bausch & Lomb, Rochester, NY 14692-0450, 1-800-553-5340.

GOOD LENS CARE PRACTICES:

- Always wash and rinse your hands before you handle your lenses.
- Clean, rinse and disinfect your lenses each time you remove them.
- Always handle the same lens, the right or the left, first in order to avoid mix-ups.
- After use, always empty and rinse the lens case with fresh rinsing solution and allow to air dry.
- When wearing disposable lenses, each time you remove a lens, discard it and replace it with a new lens.

HOW SUPPLIED:

Bausch & Lomb ReNu Preservative Free Lubricating & Rewetting Drops is supplied in 4 foil pouches each containing 5 sterile 0.01 FL. OZ. (0.3 mL) Disposable, Single-Use Containers. Cartons and container tabs are marked with a lot number and expiration date.

 Toll Free Product Information
Call: 1-800-553-5340

Bausch & Lomb and ReNu are registered trademarks of Bausch & Lomb Incorporated.

© Bausch & Lomb Incorporated

Product manufactured in Italy for:
Bausch & Lomb, Rochester, NY 14692-0450 Order No. 625224 5/99